

Contents

Preface	XI	2 Negative Casting Techniques	
Acknowledgements	XII	Introduction	17
Introduction		Neutral Position Casting Techniques	18
The Development of the 'Root' Functional Orthosis	1	Suspension Technique	18
References	4	Advantages and Disadvantages of the Suspension Technique	25
1 Indications for the Use of Functional Orthoses		Common Casting Errors When Using the Suspension Technique	25
Introduction	5	Modified Suspension Technique	26
Indications for the Use of Functional Orthoses	6	Advantages and Disadvantages of the Modified Suspension Technique	29
Support of Abnormal Structural Positions of the Forefoot	6	Common Casting Errors Using the Modified Suspension Technique	30
Forefoot Invertus	6	Direct Pressure Technique	30
Forefoot Evertus	7	Advantages and Disadvantages of the Direct Pressure Technique	30
Support of Rearfoot Deformity	8	Pronated Position Casting Techniques	32
Subtalar Joint Varus	8	Semi-Pronated Technique	32
Subtalar Joint Valgus	9	Fully Pronated Technique	33
Resistance to Abnormal Forces from the Leg	9	References	34
Extrinsic Forces Causing Abnormal Pronation of the Foot	9	Bibliography	34
External Forces Causing Abnormal Supination of the Foot	11	3 Evaluation of the Negative Cast	
Reposition of the Calcaneus	13	Introduction	35
Correcting a Medial Heel Strike	13	Bisection of the Posterior Heel of the Negative Cast	35
Correcting an Excessive Lateral Heel Strike	13	Lateral Curvature Method	36
Provision for a Normal Degree of Contact Phase Shock Absorption	14	Calcaneal Bisection Method	38
Immobilisation of the Subtalar Joint	14	Critique of a Neutral Position Cast	39
Subtalar and Midtarsal Joint Arthritis Causing Severe Pain with Joint Motion	14	Exterior Assessment of the Negative Cast Examination of the Forefoot to Rearfoot Position	44
Tarsal Coalitions	15	Interior Assessment of the Negative Cast	45
References	16	Critique of a Pronated or Semi-Pronated Cast	47
		References	48
		Bibliography	48

4 Orthosis Fabrication		
Introduction	49	
'Root Functional' versus 'Standard Biomechanical' Orthosis	50	
The Positive Cast	50	
The Orthosis Shell	51	
Forefoot and Rearfoot Control	53	
Part 1 – Negative and Positive Cast		
Preparation	54	
Preparing the Negative Cast for Pouring	54	
Preparing the Plaster for Pouring	55	
Pouring the Positive Cast	55	
Evaluation of the Forefoot to Rearfoot Position	57	
Marking the Positive Cast	58	
Marking the 1st and 5th Metatarsal Head Positions	58	
Marking the Distal Heelcup Line	60	
Marking the Heelcup Height	61	
Marking the Area for the Medial Addition	62	
Plaster Modifications to the Positive Cast	62	
Modifying the Cast for Intrinsic Forefoot Balancing	62	
Preparing the Positive Cast for Extrinsic Forefoot Posting	68	
The Application of a Medial Addition	72	
The Application of a Lateral Expansion	74	
Special Positive Cast Modifications	76	
Part 2 – Manufacture of the Orthosis Shell	80	
Preparing the Positive Cast for Shell Fabrication	80	
Marking the Outline of the Orthosis Shell on the Positive Cast	80	
Orthosis Shell Materials	81	
Pressing the Orthosis Shell	82	
Grinding the Orthosis Shell	84	
Part 3 – The Application of Forefoot and Rearfoot Posts	86	
Posting Materials	86	
Forefoot Posting the Orthosis Shell	86	
The Concept of the Forefoot Post	86	
The Dimensions of a Forefoot Post	86	
Preparing the Shell for Forefoot Posting	88	
Applying a Forefoot Post	88	
Rearfoot Posting the Orthosis Shell	91	
The Concept of the Rearfoot Post	91	
The Dimensions of a Rearfoot Post	91	
Preparing the Orthosis Shell for Rearfoot Posting	91	
Applying a Rearfoot Post	91	
'Clearing' the Posts	96	
Part 4 – Finishing the Functional Orthosis	96	
Ensuring an Accurate Heelcup Height and Orthosis Shell Dimensions	96	
Ensuring an Accurate Heelcup Height	96	
Ensuring Accurate Orthosis Shell Dimensions	97	
Smooth Grinding the Dimensions of the Forefoot Post	98	
Smooth Grinding the Dimensions of the Rearfoot Post	98	
Polishing the Orthosis	99	
Engraving the Orthosis Shell	100	
Top Covering the Orthosis Shell	100	
Part 5 – Alternative Fabrication		
Techniques	100	
Introduction	100	
Orthosis Shell Modifications	101	
Rearfoot Control	102	
Intrinsic (Cast Correction) Rearfoot Control	102	
Modified Intrinsic (Shell Modification) Rearfoot Control	104	
Polysectional Biaxial Posting (Cast Correction) for Rearfoot Control	104	
Heel Lifts	105	
Forefoot Control	106	
Polysectional Triaxial 'Posting' (Cast Correction) for Forefoot Control	106	
Computer-Aided Fabrication Processes	107	
References	107	
5 Prescription Writing		
Introduction	109	
The Biomechanical Data Required to Formulate an Accurate Orthosis Prescription	110	
Neutral Calcaneal Stance Position (NCSP)	110	

Polysectional Triaxial Forefoot Control	165
Forefoot Valgus in Association with Rearfoot Deformity	165
Plantarflexed 1st Ray	167
Compensation and Pathomechanics	167
Treatment Objectives	168
The Negative Cast	169
Root Intrinsic Forefoot Balancing Technique	169
Extrinsic Forefoot Posting Technique	169
Polysectional Triaxial 'Posting'	170
Plantarflexed 1st Ray in Association with Rearfoot Deformity	170
Rearfoot Varus	170
Compensation and Pathomechanics	170
Treatment Objectives	171
The Negative Cast	171
Treatment of Rearfoot Varus Deformities	172
Rearfoot Valgus	173
Compensation and Pathomechanics	174
Treatment Objectives	174
The Negative Cast	175
Treatment of Rearfoot Valgus Deformities	175
References	178

7 Critique of a Functional Orthosis

Introduction	179
Checking the Prescription	179
Checking the Degree of Elevation Incorporated into an Extrinsic Rearfoot Post	179
Checking the Degree of Rearfoot Posting and Pronation Skive	181
Checking the Medial and Lateral Heelcup Heights	182
Checking the Degree of Pitch Incorporated into an Extrinsic Forefoot Post	185
Checking the Degree of Extrinsic Forefoot Posting	185
Checking the Degree of Intrinsic Forefoot Control Incorporated into the Orthosis Shell	186
Checking Correct Points of Manufacture	187

8 Fitting the Functional Orthosis

Introduction	190
Pre-Therapy Advice	190
Your Orthotic Devices and How to Wear Them	190
Fitting the Orthosis to the Patient's Foot	192
The Non-Weight-Bearing Foot	192
The Weight-Bearing Foot	195
Fitting the Orthosis to the Patient's Shoe	198
Checking the Patient's Shoe	198
Checking the Orthosis to Shoe Relationship	199
Gait Analysis with the Patient Wearing the Orthosis	200

9 Common Problems with Orthosis Therapy and Their Management

Introduction	201
Shoe Accomodation Problems	201
Heel Slippage	202
Forefoot Compression	203
Orthosis Instability	204
Sagittal Plane	204
Frontal Plane	204
Transverse Plane	204
Shell Irritation	204
Front Edge Irritation	205
Lateral Heelcup Irritation or Pinching	206
Posterior Heelcup Irritation or Pinching	207
Lateral Forefoot Irritation	207
Medial Forefoot Irritation	208
Medial Longitudinal Arch Irritation	209
Plantar-Cuboid Irritation	209
Orthosis Shell Fractures	210
Orthosis Squeaking	210
Soft-Tissue Aches and Pains	210
Chronic Peroneal Muscle Affections	211
Chronic Calf Muscle Affections	211
Registered Trade Marks	212
Subject Index	213