

Table of Contents

I. Restorations on natural teeth 15

II. Conditionally removable restorations 45

III. Removable restorations 57

IV. Restorations with ZrO₂ on implants 127

Further reading 145

I. Restorations on natural teeth

Telescopic restoration with a BioHPP frame and visio.lign veneers on natural abutments 16-17

*Dr. Arne Sauer, Cospeda (DE), MDT Maxi Findeiß, Dentallabor Grüttner,
Pößneck (DE)*

Oral rehabilitation full tape Minimally invasive technique 18-19

*Peidro Puerto J., Cbo Pastor M., Haya Fernández C.,
Ribes Vallés C. Facultad de Ciencias de la Salud – Odontología – Universidad
Cardenal Herrera - CEU*

BioHPP telescopic restoration of the maxilla/mandible 20-21

*DT Mathias Göricke, OralElegance Dentallabor Horst-Dieter Kraus GmbH,
Stuttgart (DE)*

Prosthetic rehabilitation with natural BioHPP frame material 22-23

Dr. Marius Klein MSc., Pfaffenhofen (DE)

Closing gaps with a BioHPP adhesive bridge 24-25

*Dr. Sabrina Hollbach, Orthodontist
Dr. Hameister & M.Kerzendorf
MDT B.Miketta, 63654 Büdingen (DE)*

Combination restoration based on telescopes and Locator® with BioHPP 26-27

Ralph Rogmans Tandtechniek BV, Heerlen (NL)

Double-crown technology with the high-performance ceramic-strengthened polymer, BioHPP 28-29

MDT Sebastian Schuldes, MSc. (DE)

High-tech prosthesis made using the thermopress 400 30-31

Oliver Heinzmann (MDT) Heppenheim (DE)

BioHPP in the treatment of children 32-33

*Dr. Ellen Ritterbusch, Opus-DC, Ulm (DE)
MDT Florian Waschek, Kaufbeuren (DE)*

Telescopic restoration using BioHPP on natural teeth 34-35

*Dr. Klaus Schütte, Hamburg (DE)
MDT Thorsten Meyer, Praxislabor Schütte, Hamburg (DE)*

Table of Contents

I. Restorations on natural teeth

BioHPP restoration on natural teeth	36-37
<i>Dr. med. dent. Manfred Schöps, Marktoberdorf (DE)</i>	
<i>MDT Florian Waschek, Kaufbeuren (DE)</i>	
Full prosthetics – from set-up through to manufacture	38-39
<i>MDT Heinzmann, Bensheim (DE)</i>	
Manufacture of double crowns using the precision casting process	40-41
<i>Dentallabor Wolfgang Klittich, Tübingen DT Norbert Gritzmann (DE)</i>	
Maxilla case description – zirconium dioxide primary work for telescopic bridge/fully anatomical secondary HIPC restoration	42-43
<i>MDT Erich Safaric (DE)</i>	

II. Conditionally removable restorations

A new restoration using BioHPP: expansion work in the maxilla and mandible using biocompatible material	46-47
<i>Dr. Boris Peter Hrala, Zell am See (AT)</i>	
<i>MDT Jana Schmid, Traunstein (DE)</i>	
Dividing mass with BioHPP	48-49
<i>Dr. H. Hameister & Martin Kerzendorf (DE)</i>	
<i>MDT Bernd Miketta, 63654 Büdingen (DE)</i>	
Hybrid bridge on a natural tooth and an implant with BioHPP and visio.lign	50-51
<i>Dr. H. Hameister & Martin Kerzendorf</i>	
<i>Zahntechnik: MDT Bernd Miketta, 63654 Büdingen (DE)</i>	
Finding new answers to old questions	52-53
<i>MDT Carsten Krause, Calw (DE)</i>	
<i>Tagesklinik für zahnärztliche Chirurgie u. Implantologie</i>	
<i>Dr. Hackenberg, Schömberg, Prothetik Dentist Häusler, Calw-Altburg (DE)</i>	
Prosthetic reconstruction of the gingiva and teeth using a BioHPP bridge with transversal screws	54-55
<i>Burzin Khan (IND)</i>	

III. Removable restorations

Clinical Experience with a New Prosthetic Material for Implant Superstructures _____ 58-59

*Priv.-Doz. Dr. Jörg Neugebauer, Dr. Steffen Kistler,
Stephan Adler, CDT, Dr. Frank Kistler (DE)*

The prosthetic concept behind SKY fast & fixed transferred to ‚all-on-six‘ _____ 60-61

*Dr. Wilhelm Spurzem, SmilePlus, Bensheim (DE)
MDT Gunnar Dietz, Dental Design, Walldorf (DE)*

Working with bridges at the limit _____ 62-63

*MDT Markus Jedlinski, Crailsheim (DE)
Dr. Ryssel, Crailsheim (DE)*

Immediate prosthetic restoration based on implants with a fixed bridge, as per the ‚SKY fast & fixed‘ concept _____ 64-65

Dr. Ludovic Heinrich, Colmar (DE)

Polyether ether ketone (PEEK) employee on dental implant prosthesis _____ 66-67

*Cabo Pastor M, Gallud Romero L, Haya Fernández C,
Peidro Puerto J. Facultad de Ciencias de la Salud – Odontología- Universidad
Cardenal Herrera - CEU*

SKY fast & fixed restoration – transversal screws with posterior galvano copings _____ 68-69

*Dr. Armin Gräfe, Zahnärzte Wellingsbütteler Markt, Hamburg (DE)
MDT Frank Köpp Zahntechnik Wellingsbütteler Markt, Hamburg (DE)
Zirconium frame: MDT Mr. Ermoneit, Paul Tax Zahntechnik, Hamburg (DE)*

SKY fast & fixed for extensive periodontal damage to residual dentition _____ 70-71

Dentist Joachim Kanzow, Göttingen (DE)

SKY fast & fixed – for temporary and definitive restorations (BioHPP) _____ 72-73

*Clinic: Dr. Thore Kreimeyer, Hannover (DE)
Dental Technology: MDT Jürgen Tinne, Hannover (DE)*

Restoration of third-party implant systems and outdated prosthetic solutions with BioHPP and visio.lign _____ 74-75

Dentist Christian Möller, Delmenhorst (DE)

III. Removable restorations

Avoiding an external sinus lift through use of an angulated implant	76-77
--	--------------

Dr. Alois Müller, Göttingen (DE)

Restoration of an edentulous jaw with BioHPP framework and visio.lign veneering system	78-79
---	--------------

Dr. VIJAY MABRUKAR, Pune (IND)

Mr. Shashikumar Nair, Shree Sai Dental Lab, Pune (IND)

Implant-based restoration for maxilla/mandible	80-81
---	--------------

MDT. v.Schön-Angerer, Kassel (DE)

PEEK as a framework material for a 12-piece implant-supported bridge	82-83
---	--------------

Dr. med. dent. Bernd Siewert, Madrid, und DT Mario Parra, Alicante (ES)

Efficient fixed restoration following loss of the last tooth	84-85
---	--------------

Praxis am Moritzplatz Dr. Lara Müller, Augsburg (DE)

Labor Miller & Schmuck MDT Jo Miller, Augsburg (DE)

Implant-supported restoration on a high-performance polymer BioHPP frame	86-87
---	--------------

DT.Mario Parra, Alicante (ES)

and Dr. Daniel Escribano, Mazarrón (ES)

Temporary immediate restoration with customised mesostructure (crea.lign), and definitive restoration with a veneer system (visio.lign)	88-89
--	--------------

Dr. Marcus Parschau, Buchholz Buchholz in der Nordheide (DE)

DT Mario Pütz, Buchholz in der Nordheide (DE)

Restoration of an edentulous jaw with ceramic bridges on BioHPP SKY elegance abutments	90-91
---	--------------

Dr. Johannes Wolf, Zahnarzt für Oralchirurgie, Bonn (DE)

Manuel Höpfner, dental technician (DE)

Labor Kurth & Evers Bornheim (DE)

Partially removable implant-based restoration manufactured from the high-performance polymer, BioHPP	92-93
---	--------------

MDT Sebastian Schuldes, MSc., Eisenach (DE)

Restoration of an edentulous maxilla with whiteSKY zirconium implants, in accordance with the SKY fast & fixed concept	94-95
---	--------------

Dr. Stehle, Lingenfeld (DE)

III. Removable restorations

Bar-based restoration with BioHPP _____ 96-97

*Dr. Wolfram Friedrich, Greiz, MDT Anja Popp,
Labor Anja Popp Zahntechnik, Greiz (DE)*

**Implant-supported telescopic prostheses composed
of zirconium/galvano/BioHPP and zirconium/BioHPP _____ 98-99**

*Dr. Fürst, Ulrich, Attnang-Puchheim (AT)
DT John Nekola, Attnang-Puchheim (AT)*

**New restoration: zirconium oxide – abutments on implants,
secondary BioHPP with CAD/CAM milling,
cobalt chrome model casting, palate-free _____ 100-101**

MDT Sven Kirch, Dentallabor Jan Langner, Schwäbisch Gmünd (DE)

**Using the high-performance polymer, BioHPP,
for an implant-based prosthetic bar restoration _____ 102-103**

*Dentist Ingo Grübel, Herne, and MDT Martin Heinrichs, Bochum
First publication in „ZAHNTECH MAG 11/2014,
Spitta Verlag GmbH (Balingen) (DE)*

mini²SKY / Locator® – a simple method for fixing prostheses _____ 104-105

*Dr. Britta Rose, Borgentreich (DE)
Dentallabor König, Warburg (DE)*

SKY Locator® A simple approach for elder patients _____ 106-107

Dr. Vincent Roubinet, Montelimar (FR)

**Mandibular movement/function before and after mandibular
mini¹SKY overdenture by kinesiographic investigation.
Case presentation _____ 108-109**

*Dr. Alexandru Eugen David, Prof. univ. Dr. Elena Preoteasa,
Conf. univ. Dr. Catalina Murariu Magureanu*

miniSKY overdenture – microbial plaque adherence _____ 110-111

*Elena Preoteasa, Cristina Rizea, Florin Obadan
Carol Davila University of Medicine and Pharmacy, Bukarest, Rumänien*

retention.sil as silicone matrices for implant overdenture _____ 112-113

*Elena Preoteasa, Ana-Maria Cristina Tancu,
Cristina Teodora Preoteasa
Carol Davila University of Medicine and Pharmacy, Bukarest (RO)*

Simple, patient-centred restorations _____ 114-115

*Dr. Dr. Rabih Nahas, MSc, Bremen (DE)
MDT Volker Dehn, Bremen (DE)*

III. Removable restorations

TiSi.snap – the economic solution to fix overdenture	116-117
---	----------------

Dr. Janak Gandhi, Mumbai (IND)

Silicone as a matrix material – a problem - solver and a template for success?	118-119
---	----------------

MDT Erwin Klampfer -Dentale Manufaktur, Saarbrücken (DE)

mini¹SKY overdenture treatment options. Case series	120-121
---	----------------

Elena Preoteasa, Laurentiu Florica, Marina Imre (RO)

Metal-free lock-based bridge	122-123
-------------------------------------	----------------

MDT Bernhard Stühler, Nürnberg (DE)

BioHPP – removable combination-type restoration with locking pin	124-125
---	----------------

Dr. Heinz-Dieter Unger, Osnabrück (DE)

IV. Restorations with ZrO₂ on implants

Excellence in reconstruction – Zirconium vs PEEK	128-129
---	----------------

Jochen K. Alius¹, Bernd Goller²

¹Dr. med. dent., MSc., Private practice Nürnberg (DE)

²MDT, Dental laboratory, FAI GmbH, Nürnberg (DE)

Laboratory-manufactured individual abutments with a ceramic bridge restoration	130-131
---	----------------

Klaus Goldschmidt, Lingen Frau Julia Schlee,

Schlee-Dental, Bocholt (DE)

SKY fast & fixed restoration with a zirconium frame in the maxilla and a BioHPP frame in the mandible	132-133
--	----------------

Sylke Armke, Ilsede (DE)

Arash Aref, MKG am Marktplatz, Peine (DE)

MDT Christoph Dyck, dyck Dental-Technik, Lengede-Broistedt (DE)

Using bridge with ceramic veneers	134-135
--	----------------

Stephan Adler, Siegfried Weiss

Landsberg am Lech (DE)

IV. Restorations with ZrO₂ on implants

**Removable dental prosthesis manufactured
from zirconium on a bar design, secured with
three bredent snap-on and rod attachments _____ 136-137**

MDT Oliver Meyer-Buchtien (DE)

DT Wladimir Dulson (DE)

**Therapeutic splints made of BioHPP?
Yes, now also possible with CAD/CAM milling! _____ 138-139**

Dentec Solutions Network, MDT Kappert, Essen (DE)

**Metal-free telescopic restoration of dentition with
an emphasis on the organism as a whole _____ 140-141**

Dipl.-Stom. Uta Hirmer and Dr. Winfried Kather, Bad Lobenstein (DE)

**The BioHPP in the immediate loading
for the aesthetic cover of abutment on fixture _____ 142-143**

V.Giannelli, L.C.Campos, C. Baroncini, G. Ghirlanda (IT)